

# eBite

Draadloos & Oplaadbaar  
Intra-oraal lichtstelsel

Liberating dentistry, one hand at a time



Draadloos



Oplaadbaar



BiteBlock



Licht



**Algold dental**

Overveld 11G  
3848 BT Harderwijk  
tel.: +31 (0)341 42 33 85  
fax: +31 (0)341 43 22 02  
info@algold.nl  
www.algold.nl

■ Product Features



■ Specifications

Biteblock	Autoclavable Silicone material (Large and small sizes available)
LED Light	Illuminance: 5,000~8,000 lux Color temperature: 5,500K Operation life: 50,000 hours
Battery	Li-Polymer Rechargeable battery
Power Adapter	Input: 100~200VAC / 50~60Hz Output: DC 5.0V / 1.0A