


Prophylaxis


White nylon prophylaxis brushes	6
Microtufts	7
White & black bristle prophylaxis brushes	8
SMART Prophylaxis brushes	9 - 10
Standard prophylaxis polishers - with latex	11
Standard prophylaxis polishers - latex free	11
SMART Prophylaxis polishers - latex free	12 - 14
SMART Mini prophylaxis polishers - latex free	15
Prophylaxis paste - 2g cups	16
Prophylaxis paste - 50g pots	16
Sica occlusal brushes	17
Dia brushes	18
Glass ceramic burs	19
Diamond standard polishing paste	20
ALPINE Minifelts	20
PERFECTIP - Ultrasonic tips	21 - 24

White nylon prophy brushes


These shapes and profiles have been designed to meet the exact needs of the dental profession. The nylon filament is for cleaning and polishing of **natural teeth**, **gold** and **amalgam fillings**. Available as snap-on, screw type and right angle, they are ideal for removing **plaque** and **stains** and are vibration-free for patient comfort. They are designed for use with **STODDARD** Prophyl polishing paste - see page 16.




REF	JC/WNRA	JC/WNGRA	KXL/WNRA	KS/WNRA001	KS/WNRA002
Material	white nylon	white nylon extra soft	white nylon		
Ø mm	7,0	7,0	6,0	5,0	5,0
Trim mm	5,0	5,0	6,0	6,0	6,0
Shank	RA	RA	RA	RA	RA
	100	100	100	100	100


 max. 10'000 rpm






REF	JC/WNSC	JC/WNGSC	KS/WNSC001	KS/WNSC002
Material	white nylon	white nylon extra soft	white nylon	
Ø mm	7,0	7,0	5,0	5,0
Trim mm	5,0	5,0	6,0	6,0
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100	100

 max. 10'000 rpm



REF	KSO/WN001	KSO/WN002
Material	white nylon	
Ø mm	5,0	5,0
Trim mm	7,0	7,0
Shank	snap-on	snap-on
	100	100

 max. 10'000 rpm



REF	SET 3028
-----	----------

Prophy Set - 3028 contains 6 pieces of JC/WNRA and 1 piece each of KS/WNRA001 - KS/WNRA002 and MT/XWNRA
See page 116


Microtufts


Specially designed for polishing and cleaning of areas with restricted access such as **cavities, fissures and the edges of the gumline, in and around braces and implants etc.** This product can be used with or without a polishing paste.

With light pressure the bristles spread themselves and can even reach under the gumline. Microtuff brushes are designed for slow speeds to clean before bonding **fillings, inlays, onlays and crowns.** Use with a disinfectant or prophyl paste.


Dentists regularly need to clean in the smallest and confined spaces. Efficient cleaning is often difficult or even impossible with conventional brushes.


We have designed a new minibrush specifically for this task: the Microtuff! The working tuft is narrower than the shank so it is easy to reach all the previously impossible small and confined areas. The Microtuff is manufactured with nylon bristles and can be used with or without **STODDARD** Prophyl paste - see page 16.



REF	MT/XWNRA	MTS/XWNRA
Material	white nylon	
Ø mm	2,0	2,0
Trim mm	3,5	3,5
Total length mm	23,5	18,5
Shank	RA	RA
	100	100

● max. 10'000 rpm




REF	MT/XWNSC	MTS/XWNSC
Material	white nylon	
Ø mm	2,0	2,0
Trim mm	3,5	3,5
Total length	18,5	13,5
Shank mm	screw type 1,72	screw type 1,72
	100	100


● ✨ max. 10'000 rpm




White & black bristle prophy brushes


These shapes and profiles have been designed to meet the exact needs of the dental profession. The white and black bristles are for cleaning and polishing of **natural teeth**, **gold** and **amalgam fillings**. Available as snap-on, screw type and right angle, they are ideal for removing **plaque** and **stains** and are vibration-free for patient comfort. They are designed for use with **STODDARD** Prophy polishing paste - see page 16.




REF	JC/WBRA	KS/WBRA001	KS/WBRA002	KS/BBRA001	KS/BBRA002
Material	white bristles (natural hair)			black bristles (natural hair)	
Ø mm	7,0	5,0	5,0	5,0	5,0
Trim mm	5,0	6,0	6,0	6,0	6,0
Shank	RA	RA	RA	RA	RA
	100	100	100	100	100


max. 10'000 rpm



REF	JC/WBSC	KS/WBSC001	KS/WBSC002	KS/BBSC001	KS/BBSC002
Material	white bristles (natural hair)			black bristles (natural hair)	
Ø mm	7,0	5,0	5,0	5,0	5,0
Trim mm	5,0	6,0	6,0	6,0	6,0
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100	100	100

max. 10'000 rpm



REF	KSO/WB001	KSO/WB002
Material	white bristles (natural hair)	
Ø mm	5,0	5,0
Trim mm	7,0	7,0
Shank	snap-on	snap-on
	100	100

max. 10'000 rpm

SMART Prophy brushes

You can now choose the stiffness of nylon that suits your specific requirements to clean **plaque** and **stains** caused by coffee, tea, nicotine, red wine, etc. These new coloured nylon brushes are manufactured in 3 gauges of high quality flexible nylon. They are designed for use with **STODDARD** Prophy polishing paste - see page 16.


Pink = soft
Purple = medium
Blue = hard

SMART Junior cup prophy brushes



REF	JC/PKNRA	JC/PUNRA	JC/PBNRA
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	7,0	7,0	7,0
Trim mm	5,0	5,0	5,0
Shank	RA	RA	RA
	100	100	100


max. 10'000 rpm



REF SET 3026

SMART Brush Set - 3026 contains 3 pieces each of JC/PKNRA - JC/PUNRA and JC/PBNRA
See page 116


SMART Pencil prophy brushes



REF	KS/PKNRA001	KS/PUNRA001	KS/PBNRA001
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	5,0	5,0	5,0
Trim mm	6,0	6,0	6,0
Shank	RA	RA	RA
	100	100	100

max. 10'000 rpm

SMART Pencil prophy brushes



REF	KS/PKNRA002	KS/PUNRA002	KS/PBNRA002
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	5,0	5,0	5,0
Trim mm	6,0	6,0	6,0
Shank	RA	RA	RA
	100	100	100

max. 10'000 rpm

SMART Microtuff prophy brushes (with long shank)



REF	MT/PKNRA	MT/PUNRA	MT/PBNRA
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	2,0	2,0	2,0
Trim mm	3,5	3,5	3,5
Total length mm	23,5	23,5	23,5
Shank	RA	RA	RA
	100	100	100

max. 10'000 rpm

SMART Microtuff prophy brushes (with short shank)



REF	MTS/PKNRA	MTS/PUNRA	MTS/PBNRA
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	2,0	2,0	2,0
Trim mm	3,5	3,5	3,5
Total length mm	18,5	18,5	18,5
Shank	RA	RA	RA
	100	100	100

max. 10'000 rpm


SMART Junior cup prophy brushes



REF	JC/PKNSC	JC/PUNSC	JC/PBNSC
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	7,0	7,0	7,0
Trim mm	5,0	5,0	5,0
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100

max. 10'000 rpm


SMART Pencil prophy brushes



REF	KS/PKNSC001	KS/PUNSC001	KS/PBNSC001
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	5,0	5,0	5,0
Trim mm	6,0	6,0	6,0
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100

max. 10'000 rpm

SMART Pencil prophy brushes



REF	KS/PKNSC002	KS/PUNSC002	KS/PBNSC002
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	5,0	5,0	5,0
Trim mm	6,0	6,0	6,0
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100

max. 10'000 rpm

SMART Microtuft prophy brushes (with long shank)



REF	MT/PKNSC	MT/PUNSC	MT/PBNSC
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	2,0	2,0	2,0
Trim mm	3,5	3,5	3,5
Total length mm	18,5	18,5	18,5
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100

max. 10'000 rpm

SMART Microtuft prophy brushes (with short shank)




REF	MTS/PKNSC	MTS/PUNSC	MTS/PBNSC
Material & colour	pink nylon	purple nylon	dark blue nylon
Stiffness	soft	medium	hard
Ø mm	2,0	2,0	2,0
Trim mm	3,5	3,5	3,5
Total length mm	13,5	13,5	13,5
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100

max. 10'000 rpm

Standard prophy polishers - with latex

These shapes and profiles have been designed to meet the exact needs of the dental profession. The rubber is for cleaning and polishing of **natural teeth, gold** and **amalgam fillings**. Available as snap-on, screw type and right angle, they are ideal for removing **plaque** and **stains** and are vibration-free for patient comfort. They are designed for use with **STODDARD** Prophyl polishing paste - see page 16.




REF	TPR003	TPY003A	TPY002	TPR003A	TPR002A	TPR004A
Material	natural rubber					
Stiffness	medium	medium	medium	medium	medium	medium
Shank	RA	RA	RA	RA	screw 1,72	snap-on
	100	100	100	100	100	100
	●	●	●		● ✦	● ✦

max. 10'000 rpm

Standard prophy polishers - latex free

The original standard prophyl polishers are now available latex free. The special medical grade syntetic polymer is ideal for cleaning and polishing of **natural teeth, gold** and **amalgam fillings**. Available as snap-on, screw type and right angle, they are perfect for removing **plaque** and **stains** and are vibration-free for patient comfort. Designed for use with **STODDARD** Prophyl polishing paste - see page 16.

- ◆ New improved
- ◆ Now available latex free



REF	EPR003	EPY003A	EPY002	EPR003A	EPR002A	EPR004A
Material	medical grade syntetic polymer - latex free					
Stiffness	medium	medium	medium	medium	medium	medium
Shank	RA	RA	RA	RA	screw 1,72	snap-on
	100	100	100	100	100	100
	●	●	●		● ✦	● ✦

max. 10'000 rpm

- Also available in the SMART range - see pages 13 - 14.
- ✦ Adapters and mandrels - see page 112.

SMART Prophy polishers - latex free



You can now choose the stiffness of the polisher that best suits your specific requirements to clean stains caused by coffee, tea, nicotine, red wine, etc. These new coloured prophy polishers are manufactured in 3 types of medical grade synthetic polymer. They are designed for use with **STODDARD** prophy polishing paste - see page 16.

Latex allergy is a medical term for a range of allergic reactions to proteins found in natural rubber latex. Due to the widespread use of latex products in healthcare and medical devices, latex allergy presents a significant healthcare problem. It is estimated that approximately 8% of the global population is sensitive to latex. The risk of a potentially life-threatening allergic reaction by a patient in the surgery has been deemed to be sufficiently high to replace all latex surgical products with synthetic alternatives.



SMART Prophy polishers - latex free (6 web)

You can now choose the stiffness of polisher that best suits your specific requirements to clean **plaque** and **stains** caused by coffee, tea, nicotine, red wine, etc. These new coloured prophy polishers are manufactured in 3 types of medical grade synthetic polymer. They are designed for use with **STODDARD** Prophy polishing paste - see page 16.


Pink = soft
Purple = medium
Blue = hard



REF	PRA001	PRA002	PRA003
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank	RA	RA	RA
	100	100	100

max. 10'000 rpm

New



REF	PRA101	PRA102	PRA103
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank	medical grade plastic shank - RA		
	100	100	100

max. 10'000 rpm

- ◆ New improved
- ◆ High quality medical plastic shank
- ◆ Latex free
- ◆ Highly economical
- ◆ High quality



REF	PSRA001	PSRA002	PSRA003
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank	RA	RA	RA
	100	100	100

max. 10'000 rpm

New



REF	SET 3027
-----	----------

SMART Polisher Set - 3027 contains 3 pieces each of PSRA001 - PSRA002 and PSRA003
See page 116



REF	PSC001	PSC002	PSC003
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100

max. 10'000 rpm



REF	PSO001	PSO002	PSO003
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank	snap-on	snap-on	snap-on
	100	100	100

max. 10'000 rpm

◆ Adapters and mandrels - see page 112.

SMART Prophy polishers - latex free (4 web and ribbed)

You can now choose the stiffness of polisher that best suits your specific requirements to clean **plaque** and **stains** caused by coffee, tea, nicotine, red wine, etc. These new coloured prophy polishers are manufactured in 2 types of medical grade synthetic polymer. They are designed for use with **STODDARD** Prophy polishing paste - see page 16.

Pink = soft
Purple = medium
Blue = hard



REF	PSRA011	PSRA012	PSRA013
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank	RA	RA	RA
	100	100	100

max. 10'000 rpm



REF	PSC011	PSC012	PSC013
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100

max. 10'000 rpm



REF	PSO011	PSO012	PSO013
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank	snap-on	snap-on	snap-on
	100	100	100

max. 10'000 rpm



SMART Mini prophy polishers - latex free (4 web and ribbed)

The SMART MINI Prophy polishers are a shorter version of our latex free products. These cups provide fast and efficient polishing where the space is restricted, specially designed to clean and polish **children's teeth** or **any other restricted areas**. These products are approximately half the size of our standard range. They are designed for use with **STODDARD** Prophylaxis polishing paste - see page 16.


Pink = soft
Purple = medium
Blue = hard



REF	PSMRA001	PSMRA002	PSMRA003
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank	RA	RA	RA
	100	100	100

max. 10'000 rpm


Comparisons between SMART Standard prophy polishers RA and SMART Mini prophy polishers RA:



REF	PSMSC001	PSMSC002	PSMSC003
Material	medical grade synthetic polymer		
Colour	pink	purple	blue
Stiffness	soft	medium	hard
Shank mm	screw type 1,72	screw type 1,72	screw type 1,72
	100	100	100

max. 10'000 rpm

Comparisons between SMART Standard prophy polishers screw type and SMART Mini prophy polishers screw type:



REF SET 3029

SMART Mini Prophylaxis Set - 3029 contains 3 pieces each of PSMRA001, PSMRA002 and PSMRA003
See page 116

Prophylaxis

Stoddard

Prophylaxis paste - 2g cups

STODDARD Prophylaxis paste in easy 2g cups is available in a medium and coarse blend of mint flavoured polishing paste formulated to have an easy to use consistency without causing splatter. Contains 1.23 % fluoride.



- ◆ Non-splatter formula
- ◆ Contains 1.23 % fluoride
- ◆ Mint flavoured
- ◆ Easy to rinse
- ◆ Gluten free
- ◆ Available in 2 grits

**Contains
1.23 %
fluoride**

1 box contains

- ◆ 200 unit cups
- ◆ 1 autoclavable prophy ring



REF	PAS007	PAS008	PAS/FR001
Flavour	mint		Autoclavable prophy ring .
RDA Value	170	220	
Grit & colour	medium - green	coarse - green	
	200 cups and 1 prophy ring	200 cups and 1 prophy ring	10 pcs.

Prophylaxis paste - 50g pots

STODDARD Prophylaxis paste is available in a medium blend of mint flavoured polishing paste formulated to have an easy to use consistency without causing splatter. This product is made of natural ingredients and does not contain fluoride. Available in an economical 50g pot.



REF	PAS001
Flavour	mint
RDA Value	170
Grit & colour	medium - pink
	1 pot 50g


**100 %
natural
ingredients**

Sica occlusal brushes

STODDARD Sica occlusal polishing brushes provide a single-step, fast and effective polishing system for **restorations** and **natural teeth**. Polishing paste is not necessary as the special filaments are loaded with fine silicon carbide polishing particles for an exceptionally high lustre finish. As the brushes wear, abrasive grit is consistently exposed on each filament at the polishing surface.

High quality results are achieved more efficiently than with conventional polishers. **STODDARD** Sica occlusal brushes are available in three shapes for easy polishing of occlusal and interproximal areas: cup, point and pencil.

- ◆ A high lustre satin finish polishing of restorations made of composites, compomers and glass ionomer cements
- ◆ No polishing paste required
- ◆ One step finishing and polishing
- ◆ Versatile
- ◆ Faster and more efficient than polishing paste systems
- ◆ Very economical
- ◆ Available with gold cups
- ◆ Low heat generation
- ◆ Autoclavable
- ◆ Reuseable



REF	JC/SICRA/B	KS/SICRA/C	KS/SICRA/D	JC/SICRA/G	KS/SICRA/E	KS/SICRA/F
Material	silicon carbide impregnated nylon					
Ø mm	7,0	5,0	5,0	7,0	5,0	5,0
Trim mm	5,0	6,0	6,0	5,0	6,0	6,0
Shank	RA	RA	RA	RA	RA	RA
Type	standard			with gold plated cup		
	10	10	10	10	10	10

max. 5'000 rpm

This is a multi-use product!



REF SET 3001

Occlusal Brush Set I - 3001 contains 3 pieces each of JC/SICRA/B - KS/SICRA/C and KS/SICRA/D
See page 117



REF SET 3004

Occlusal Brush Set II - 3004 contains 3 pieces each of JC/SICRA/G - KS/SICRA/E and KS/SICRA/F
See page 117

● Also available with HP shank in the laboratory section - see page 72.

Prophylaxis

Dia brushes



Stoddard

STODDARD Dia brush is a single step highly versatile fine polishing system for all kinds of restorative materials. Polishing paste is not required as the filaments are loaded with fine diamond polishing particles for an exceptionally high lustre finish. Diamond grit is consistently exposed on each filament at the polishing surface.


High quality results are achieved more efficiently than with conventional polishers. **STODDARD** Dia brushes are available in three shapes for easy polishing of occlusal and interproximal areas: cup, point and pencil.

- ◆ A high lustre satin finish polishing
- ◆ No polishing paste required
- ◆ One step finishing and polishing
- ◆ Versatile
- ◆ Faster and more efficient than polishing paste systems
- ◆ Very economical
- ◆ Low heat generation
- ◆ Autoclavable
- ◆ Reuseable

- ◆ All composites
 - Macrogrits
 - Hybridgrits
 - Nanogrits
 - Glasionomere
 - Compomer

- ◆ All ceramics
 - CEREC
 - Zirconia
 - All ceramic

- ◆ Amalgam
- ◆ All precious and non-precious metals
- ◆ Natural teeth



REF	JC/DIA/A	KS/DIA/A/001	KS/DIA/A/002
Material	diamond impregnated nylon		
Ø mm	7,0	5,0	5,0
Trim mm	5,0	6,0	6,0
Shank	RA	RA	RA
	10	10	10

max. 5'000 rpm

This is a multi-use product!

REF SET 3022

Dia Brush Set - 3022 contains 3 pieces each of JC/DIA/A - KS/DIA/A/001 and KS/DIA/A/002
See page 117





● Also available with HP shank in the laboratory section - see page 72.

Glass ceramic burs

STODDARD Multipurpose glass ceramic burs can be used to remove **cement, composites, adhesives** and **stains** without risk to the enamel, ceramic, metal and soft tissues. The glass ceramic bur is designed for the removal of excess only - for the initial removal of thicker layers we recommend using our carbide bur C23R/RA/016 on page 33 followed by the glass ceramic burs.

They can also be used in **orthodontics, periodontology, implantology** and **clinical process procedures in general**. Available in 6 different shapes.



REF	FB001	FB002	FB003	FB004	FB005	FB006
Material	fibreglass					
L mm	11,0	11,0	11,0	10,0	6,0	5,0
Ø mm	2,3	3,3	1,8	2,3	2,2	3,3
Shank	RA	RA	RA	RA	RA	RA
	5	5	5	5	5	5

max. 10'000 rpm



REF SET 1015


Glass ceramic Bur Set - 1015 contains 1 piece each of FB001 - FB002 - FB003 - FB004 - FB005 and FB006
See page 117

Prophylaxis

Stoddard

Diamond standard polishing paste




REF	PASD002
Content	2,0 g
Grit	1 - 2 µm
	1 syringe

The diamond standard polishing paste, 2,0g extra fine diamond grit, 1-2 microns produces an excellent high gloss and natural brightness on **all restoration materials, composites, ceramics, gold etc. and tooth enamel.**

The grey coloured (flavourless) paste is specially formulated for intra and extra-oral use, hydro-soluble for easy cleaning.

ALPINE Minifelts



REF	FP001	FP002	FP003
Material	hard felt		soft felt
L mm	3,0	5,0	3,5
Ø mm	9,0	9,0	12,0
Shank	RA	RA	unmounted
	50	50	50

The ALPINE Minifelt, the disposable felt polisher for an exceptionally high shine result (single use). Designed for use with **STODDARD** Diamond polishing paste.

Finishing/reglazing takes 30 seconds or less.

✦
↻ max. 5'000 rpm